

# AVENIO ctDNA and Tumor Tissue Targeted Panel



This 17-gene panel, which comes with the AVENIO ctDNA Targeted Kit and the AVENIO Tumor Tissue Targeted Kit, contains biomarkers referenced in guidelines, including those published by the U.S. National Comprehensive Cancer Network (NCCN).\*

| Gene                 | All Coding Regions | CRC            | Lung | Breast | Gastric | Glioma | Melanoma | Ovarian | Thyroid | Pancreatic | Other Cancers  | SNV | Indel <sup>‡</sup> | Fusion <sup>§</sup> | CNV <sup>§</sup> |
|----------------------|--------------------|----------------|------|--------|---------|--------|----------|---------|---------|------------|--|-----|--------------------|---------------------|------------------|
| ALK                  |                    |                | •    |        |         |        |          |         |         |            | ALCL, Large B-Cell Lymphoma, Neuroblastoma, Inflammatory Myofibroblastic Tumor, Spitzoid Tumor, Rhabdomyosarcoma     | •   | •                  | •                   |                  |
| APC                  |                    | •              |      |        | •       | •      |          |         |         | •          | HCC, Medulloblastoma, Desmoid  | •   | •                  |                     |                  |
| BRAF                 |                    | •              | •    |        |         | •      | •        | •       | •       |            | Borderline Ovarian, Cholangiocarcinoma, Spitzoid Tumor, Pancreas Acinar Carcinoma, Melanocytic Nevus, GIST           | •   | •                  |                     |                  |
| BRCA1                | •                  |                |      | •      |         |        |          | •       |         | •          |  | •   |                    |                     |                  |
| BRCA2                | •                  |                |      | •      |         | •      |          | •       |         | •          |  | •   |                    |                     |                  |
| DPYD**               |                    | •              |      |        |         |        |          |         |         |            |  | •   |                    |                     |                  |
| EGFR                 | •                  |                | •    |        |         | •      |          |         |         |            |  | •   | •                  |                     | •                |
| ERBB2                | •                  | • <sup>†</sup> | •    | •      | •       | •      |          | •       |         |            |  | •   | •                  |                     | •                |
| KIT                  |                    |                |      |        |         |        | •        |         |         |            | AML, Testicular Germ Cell Tumors, GIST, Thymic Carcinoma   | •   | •                  |                     |                  |
| KRAS                 | •                  | •              | •    |        | •       |        |          | •       | •       | •          | AML, Juvenile Myelomonocytic Leukemia  | •   |                    |                     |                  |
| MET                  | •                  |                | •    |        |         |        |          |         |         |            | Head and Neck Squamous Cell, HCC, Papillary RCC  | •   | •                  |                     | •                |
| NRAS                 |                    | •              | •    |        |         |        | •        |         | •       |            | MM, AML, Juvenile Myelomonocytic Leukemia  | •   |                    |                     |                  |
| PDGFRA               |                    |                |      |        |         |        |          |         |         |            | GIST, Hypereosinophilic Syndrome   | •   |                    |                     |                  |
| RET                  |                    | •              | •    |        |         |        |          |         | •       |            | Pheochromocytoma, Spitzoid Tumor, Mucosal Neuroma  | •   |                    | •                   |                  |
| ROS1                 |                    |                | •    |        | •       | •      |          |         |         |            | Spitzoid Tumor, Cholangiocarcinoma, Borderline Ovarian   | •   |                    | •                   |                  |
| TP53                 | •                  | •              | •    | •      | •       | •      |          | •       |         | •          | Various Leukemias, Adenocortical, Basal Cell Carcinoma, Esophageal, Soft Tissue and Bone Sarcomas, Head and Neck SCC | •   | • <sup>§§</sup>    |                     |                  |
| UGT1A1 <sup>††</sup> |                    | •              |      |        |         |        |          |         |         |            |  | •   |                    |                     |                  |

All coding regions are based on the longest transcript from Ensembl build 82.  
 \* National Comprehensive Cancer Network. <http://www.nccn.org>. October 15, 2016.  
 \*\*Inactivating mutations are associated with sensitivity to fluoropyrimidines.  
 † Kavuri et al. Cancer Discovery 8: 832-41 (2015). PMID: 26243863.  
 †† Certain variants are associated with sensitivity to irinotecan.  
 ‡ Indels are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest", except for EGFR exon 19 long deletions, EGFR exon 20 long insertions and MET long insertions.  
 § Detection of Fusions and CNVs are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest" in the AVENIO analysis software.  
 §§ Indel for TP53 is called for tissue only, not plasma.

台灣羅氏醫療診斷設備股份有限公司  
 104台北市中山區民權東路三段2號10樓  
 (02) 2183 6688  
 MC-TW-01068  
 僅供研究用 · 不可用於臨床診斷

AVENIO is a trademark of Roche.  
 © 2018 Roche Sequencing Solutions, Inc.  
 SEQ100325 09/2018

Associations of variations or mutations in these genes with the listed diseases were sourced from the following: COSMIC, GeneCards, My Cancer Genome. Analysis performed 11/2016.

**For Research Use Only. Not for use in diagnostic procedures.**